PaX-i3D (PHT-60CF0)

User Manual | Version 1.3.0

ENGLISH

FRANÇAIS

PORTUGUÊS

ESPAÑOL

For the American (North, South) Region

ENGLISH

Notice	05
1 PaX-i3D Green Imaging System Overview	07
1.1 General View of the PaX-i3D Green	08
1.2 Control Panel	10
1.3 Touch Screen.....	11
2 Getting Started	15
2.1 Turning on the PaX-i3D Green	15
2.2 Running the Image Viewer.....	16
2.3 Initiating the Imaging Program	20
3 Acquiring PANO Images	22
3.1 Setting Exposure Parameters	22
3.2 Positioning the Patient	26
3.3 Initiating X-Ray Exposure	38
4 Acquiring CEPH Images	40
4.1 Setting the Exposure Parameters	40
4.2 Positioning the Patient	44
4.3 Initiating X-Ray Exposure	51
5 Acquiring CT Images.....	53
5.1 Setting the Exposure Parameters	53
5.2 Positioning the Patient	58
5.3 Acquiring a SCOUT Image	63
5.4 Initiating X-Ray Exposure	66

Notice

Thank you for purchasing the **PaX-i3D Green** extra-oral imaging system.

The user manual is a part of product.

This manual describes how to operate the **PaX-i3D Green** system. It is recommended that you thoroughly familiarize yourself with this manual in order to make the most effective use of this equipment. Observe all cautions, safety messages and warnings which appear in this manual.

The **QR** code linked to video file about image capture for the **PANO Standard** and **CEPH Lateral** is provided in the manual. The smart phone and pad which have the QR code reader application program can be used to watch video.

The information contained in this manual is subject to change without prior notice to the persons concerned. For the most current information, please contact us at:

VATECH Co., Ltd.

Phone: +82-1588-9510

E-mail: gcs@vatech.co.kr

Website: www.vatech.co.kr

Manual Name: PaX-i3D Green (Model: PHT- 60CFO) User Manual

Version: 1.3.0

Publication Date: 2014-4

Copyright © 2014 by VATECH

All rights reserved.

The documentation, brand name and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations which may be required due to technical improvement.

Conventions in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions which accompany it.

To prevent personal injury and/or damage to the equipment, please observe all warnings and safety information included in this document.

	WARNING	Indicates information that should be followed with the utmost care. Failure to comply with a warning may result in severe damage to the equipment or physical injury to the operator and/or patient.
	CAUTION	Indicates a situation that demands prompt and careful action, a specific remedy, or emergency attention.
	X-RAY	Indicates a possible danger of exposure to radiation.
	IMPORTANT	Indicates a situation or action that could potentially cause problems to the equipment and/or its operation.
	NOTE	Emphasizes important information or provide useful tips and hints.
	SINGLE USE	Indicates a component which must be replaced for each new patient.

1 PaX-i3D Green Imaging System Overview

The **PaX-i3D Green** is an advanced digital dental diagnostic system that incorporates Panoramic (PANO), Cephalometric (CEPH), and 3D CT imaging capabilities into a single system.

The **PaX-i3D Green** system is intended for use to take a panoramic, cephalometric and CT image of the oral and craniofacial anatomy to provide diagnostic information for adult and pediatric patients.

The **PaX-i3D Green** may only be used by dentists, X-Ray technicians and other professionals who are licensed to perform X-Rays by the law of the region in which it is being used.

Standards and Regulations

The PaX-i3D Green is designed and manufactured to meet the following standards:

IEC/EN/UL 60601-1, IEC/EN 60601-1-1, IEC/EN 60601-1-2, IEC/EN 60601-1-3, IEC/EN 60601-2-7, IEC/EN 60601-2-28, IEC/EN 60601-2-32, IEC/EN 60601-2-44, ISO 9001, ISO 13485

NEMA standard publication PS 3.1-3.18, 2008

X-Ray source assembly [HDG-07B10T2] IEC 60601-2-28 (1993)

The CE symbol grants this equipment compliance with the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.

UL symbol grants this equipment compliance with the UL 60601-1 / CAN/CSA C22.2 No.601.1


Classifications (IEC60601-1 6.1):

Protection against the ingress of water: Ordinary Equipment: IPX0

Protection against electric shock: Class 1 equipment, Type B Applied Parts

1.1 General View of the PaX-i3D Green


No.	Item	Details
1	X-RAY SENSOR (CEPH)	CEPH imaging sensor (optional). - 1-1. One shot Type - 1-2. Scan Type
2	NASAL POSITIONER	- Positions the patient during CEPH imaging. - The ruler for use to compensate the acquired image that is different from actual size.
3	EAR ROD	Secures the patient's head during CEPH imaging.
4	SECONDARY COLLIMATOR	Limits the X-Ray irradiation field for CEPH scanning.
5	X-RAY TUBE	The source of X-Ray emission.
6	TEMPLE SUPPORT	Supports the patient's temples during PANO scan.
7	ROTATING UNIT	Rotates around the patient's head during exposure (depending on the imaging mode).
8	X-RAY SENSOR (PANO / CT)	Sensors: used for CT and PANO imaging. The sensor will be switched automatically (CT ↔ PANO, and vice versa).
9	CHINREST	Supports the patient's chin during image acquisition.
10	CONTROL PANEL	Allows the operator to control certain functions of the equipment and displays operational parameters. For more details, refer to 1.2 Control Panel and 1.3 Touch Screen
11	HANDLE FRAME	Held firmly by the patient during imaging to stabilize his / her position.
12	TELESCOPIC COLUMN	Moves the column up or down for patient positioning.
13	STATIONARY COLUMN	The portion of the column fixed to the base unit.
14	BASE	Used to balance and stabilize the equipment.

NOTE

The structural overview of the PaX-i3D Green may differ depending on the model.

1.2 Control Panel

Buttons		Functions
1	Touch Screen	Configures the parameter settings in each imaging mode. For more information on this, refer to 1.3 Touch Screen .
2	Laser Beam ON / OFF	Turns ON or OFF the laser beams.
3	READY	Indicates that parameter settings and the patient alignment are finished and imaging is ready.
4	Chinrest UP	Moves the chinrest up. - Functions only in the CT imaging mode.
	Frankfurt plane laser beam UP	Moves Frankfurt plane laser beam up. - Functions only in the PANO imaging mode.
5	Column UP	Moves the column up.
6	Chinrest DOWN	Lowers the chinrest. - Functions only in the CT imaging mode.
	Frankfurt plane laser beam DOWN	Moves Frankfurt plane laser beam down. - Functions only in the PANO imaging mode.
7	Column DOWN	Lowers the column.

1.3 Touch Screen

You can configure the parameter settings in each mode using the touch screen. The **touch screen** on the equipment and **the imaging program** on the PC are synchronized in real time, and display the same parameter settings.


Never use pointed objects, such as a ballpoint pen or pencil, to manipulate the LCD screen as they could scratch or damage the LCD screen.

CBCT


No.	Function	Description
1	Imaging parameter settings information	Displays currently set information such as patient type and X-Ray setting values.
2	FOV selection	Shows the CT sensor model mounted.
3	Selection settings	Displays the current selections in sequential order.
4	Imaging mode selection	Select the imaging mode.
5	INFO	Enables the user to set a patient's gender and X-ray intensity and controls kVp / mA.
6	Back	Moves to the previous level
7	Forward	Moves to the next level
8	CONFIRM	Click this button when the parameter settings are finished.


No.	Function	Description
1	Examination area selection	Selects the tooth position
2	Guidance Image	—

INFO


No.	Function
1	Tube voltage UP / DOWN
2	Tube current UP / DOWN
3	Patient's gender selection
4	Patient's X-ray selection
5	Displays the parameter settings

Temple Support Adjustment Button


RETURN

Returns the Rotating Unit back to its initial position.

PANO

No.	Function	Description	
1	Imaging parameter settings information	—	
2	Pano Type	Magic Pan	Optional
3	Image Quality	UHD	Optional
4	Arch Selection	Selects the type of patient's dental arch	
5	Orthogonal mode	Minimizes overlapping in the image	
6	Selection settings	Displays the current selections in sequential order.	

CEPH: OP Model (Oneshot type)

NOTE

Exposure time can be adjusted by resolution of 0.1 s in the range of 0.5 s to 1.2 s (One shot type only).

CEPH: SC Model (Scan Type)

No.	Function	Description	
1	Imaging parameter settings information	—	
2	CEPH Examination	—	
3	CEPH FOV Size (inches)	—	
4	Exposure Time Adjustment	OS / OP models	Enable the user to adjust the exposure time with the UP and DOWN buttons.
5	Image Quality	Normal	Optional
		Fast	Fast Scan

2 Getting Started

2.1 Turning on the PaX-i3D Green




Do not position the patient in the unit while it is initiating. The patient could be injured if the equipment malfunctions.

- Condensation could form inside the equipment if it is a different temperature than the surrounding room. Only turn on the equipment once it has reached room temperature.
- Wait at least 20 seconds after the equipment has been turned off to turn it back on.
- Allow the equipment to warm up for at least 5 minutes before acquiring images or preferably more than 30 minutes for image quality.

1. Turn the **ON / OFF** switch underneath the handle frame to the **ON** position.

2. Insert the license key for **Ez3D plus** into the USB port on the back of PC. The license key for **Ez3D plus** must be attached to view or analyze the 3D images.

2.2 Running the Image Viewer

EasyDent / EzDent-i is a basic imaging platform for all VATECH's dental X-Ray equipments. The Imaging Program is interfaced with **EasyDent / EzDent-i**. On your desktop, double-click **EasyDent / EzDent-i** icon. The **EasyDent / EzDent-i** main window will be displayed.

NOTE


For further details on this subject, refer to the **EasyDent / EzDent-i** user manual.

2.2.1 Creating a New Patient Record

To create a new patient record, follow the procedure outline below:

EasyDent

1. Click the **Patient** icon () on the upper left corner of the EasyDent's main GUI window.

2. The following dialog box will open.

3. Enter the required patient information. **Chart Number, First Name, and Last Name** are required fields which must be filled in. All other fields are optional, but it is recommended that they be filled in.

4. Click **Add** to save the patient record.

EzDent-i

1. Click the **Add Patient** icon from main GUI window.

2. Enter the required patient information. The **Chart Number, E-Mail address, First Name, and Last Name** are required fields which must be filled in. (The chart number fills in automatically.)


3. Click **Add** to save the patient record.

2.2.2 Retrieving Patient Records

EasyDent


You can search through the patient database using a patient's chart number, their first name, or their last name.

1. On the **Patient information pane**, double-click the **Chart No.**, **First name**, or **Last name** of the patient and the virtual keyboard will pop up.


2. Enter the **Chart No.**, **First name**, or **Last name** of the patient by clicking the mouse on the virtual key board and click the **Enter**.

3. Patient information can be displayed on the **Patient information pane** and **Patient List**.

EzDent-i


1. Enter the name or chart number of the patient to be searched on the **Search** pane and then click the **Search** button. The information on the patient that fits the search condition appears.

Double-click the Keyboard icon to display the virtual keyboard. You may search patient information using the virtual keyboard.


2. Double-click the patient information to see more details about the patient as shown below.


2.3 Initiating the Imaging Program

EasyDent

1. First, click the patient information in the patient list, and click the **Dental CT** icon () in the upper left corner of the EasyDent's main window to open the imaging program.


2. The following imaging program window opens. The sole purpose of this window is to control equipment settings and acquire images.


EzDent-i

1. Search and select the patient to be captured.

2. Click **Acquisition** and the imaging mode (CT, Panorama, or Cephalo).

3. The main screen for the selected mode appears. From the main screen, you can configure the imaging parameter settings prior to acquiring an image.

Please proceed to the next chapter.

NOTE

Refer to chapters 3 - 5 for information regarding image acquisition.

3 Acquiring PANO Images

3.1 Setting Exposure Parameters

Perform the following procedures to select the capture parameters for the specific patient and capture mode.

NOTE

You can set the imaging parameters on either the touch screen or the imaging program running on the PC. They are synched in real time and display the same environmental settings.

1. Choose an imaging mode under **Imaging Mode**.


2. Select the PANO Type.

Mode	Details	
Normal	Normal image.	
Magic PAN	Image with ultra-high resolution	optional

3. Select the image quality.

Mode	Details
HD	Image with higher resolution than Normal . Takes longer scan time than the Normal .
Normal	Normal image


4. Select the patient's arch type. By default, **Normal** is selected.

Orthogonal: This mode enables overlapping regions of teeth to be minimized when acquiring images in the ROI. If orthogonal arch is selected, its sub-modes are activated.

5. Select the ROI for panoramic image acquisition under **PANO Examination**.

6. The patient's gender will be automatically selected according to the patient information registered in EasyDent / EzDent-i. Nevertheless, ensure that this information is correct.

Age Group		VATECH's Standard
Adult	Child	≤ 12
	Man	≥ 13
	Woman	

7. Select X-ray intensity.


NOTE

X-ray intensity (Hard, Normal, Soft) is according to the operator's decision.

Soft \leq Normal \leq Hard

8. A default value for the tube voltage (kVp) and current (mA) will be displayed based on the gender and X-ray intensity. If necessary, further adjustments can be made using the arrows to the right of each number.

9. Click the **CONFIRM** button for these parameters to take effect.

NOTE

Please wait for a moment, while the rotating unit moves to its initial scanning position.

When you click the **CONFIRM** button,

- The **READY** button will begin blinking to show that it has been activated.
- The rotating unit will move to its initial scanning position.
- Three laser beams, **Mid-sagittal plane**, **Frankfurt plane** and **Canine**, will be activated to facilitate patient positioning. These beams will disappear after 20 minutes or if the **READY** button is clicked.

NOTE

Laser beam On / Off buttons:

: Imaging Program

: Control Panel

Scan Time : 0.0
DAP : 0.00 uGy x m^2


- The scan time and DAP (Dose Area Product) values will be shown on the main display for the selected exposure.

10. Guide the patient to the equipment. Position the patient within the equipment. For further information about patient positioning, refer to section **3.2 Positioning the Patient**.

3.2 Positioning the Patient

Before Positioning the Patient

- Ask the patient to remove all jewelry and metallic objects, such as earrings, hair pins, spectacles, dentures, and orthodontic appliances. These items may cause shadow images which could obscure diagnosis.
- It is strongly recommended that the patient wear a lead apron for protection against any possible scatter radiation.
- Adjust the column height to the patient's height using the up/down button or switch.

Correct positioning is an important factor in capturing the best possible image. Proper positioning will reduce the appearance of the cervical spine in the image.

Disinfect the chin support and bite block using an alcohol-based cleaning solution and wipe away all residue with dry cloth.

If you don't adjust the laser beams for patient alignment (vertical, horizontal, and canine beams), part of the image may be zoomed in or out or a ghost may appear in the image.


Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.

3.2.1 PANO Standard

Positioning the patient

*To link to the video-guided manual: scan the QR code with smart phone or smart pad.

1. Insert the bite block into the Chinrest receptacle on the equipment.
2. Place a hygienic cover over the bite block.


WARNING

The hygiene cover for the bite block and chin support is intended for single use only. Always replace the hygiene cover for each new patient.


CAUTION

Disinfect the Chinrest, Chin support and bite block using an alcohol-based cleaning solution and wipe away all residue with dry cloth.

3. Open the temple supports by clicking the temple support adjustment button.


4. Guide the patient —facing the Chinrest — to the equipment.
5. Adjust the height of the column using the column up/down button or switch (optional) until the patient's chin rests comfortably on the Chinrest.

6. Ask the patient to:

- Stand upright
- Firmly grip the handle
- Lean his/her chest lightly against the equipment
- Position his/her feet slightly forward and open stance

7. Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright.

8. Have the patient bite the bite block along its grooves using his/her upper and lower incisors

9. Secure the bite block firmly by turning the locking knob to the right.

10. Have the patient

- Close his/her lips around the bite block
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes


Ask the patient to remain still until scanning is completed.


To acquire the best image possible, ask the patient not to:

- Breathe or swallow saliva during image acquisition
- Move during image acquisition

Positioning the edentulous patient

1. Remove the bite block from the Chinrest receptacle on the equipment.


2. Open the temple supports by clicking the temple support adjustment button.
3. Guide the patient —facing the Chinrest— to the equipment.

4. Adjust the height of the column using the column up/down button or switch (optional) until the patient's chin rests comfortably on the Chinrest.

5. Ask the patient to:

- Stand upright
- Firmly grip the handle
- Lean his/her chest lightly against the equipment
- Position his/her feet slightly forward and open stance

6. Have the patient

- Close his/her lips around the bite block
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes


Ask the patient to remain still until scanning is completed.

Aligning the Laser Beams


Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.



- 1. Mid-sagittal plane laser beam:** Position the Mid-sagittal plane laser beam on the center of the patient's face to prevent magnifications of the left or right side in the final image.
- 2. Frankfurt plane laser beam:** Position the patient's head so that the Frankfurt plane is aligned with the Frankfurt plane laser beam.

- 3. Canine laser beam:** Have the patient smile to properly position the canine laser beam on the center of the patient's canine tooth.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head by clicking the temple support adjustment button.

2. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

3. Proceed to section **3.3 Initiating X-Ray Exposure**.

3.2.2 TMJ Open


There are two TMJ sub modes: **TMJ Open** and **TMJ Close**. Scan a TMJ Open image first and then scan a TMJ Close image. To prepare and position the patient, follow the process outlined below.

For a child with a large head circumference, Man/Woman mode may be selected instead of Child.

Positioning the patient


Chin support
Lock / release knob
for bite block

1. Insert the Chin support into the Chinrest receptacle on the equipment.


Disinfect the Chinrest using an alcohol-based cleaning solution and wipe away all residues with dry cloth before proceeding any further.

2. Open the **Temple Supports** by clicking the temple support adjustment button.

3. Guide the patient —facing the Chinrest — to the equipment.
4. Adjust the height of the equipment using the column up/down button or switch (optional) until the patient's chin is resting on the Chinrest. Ensure that the chin is in contact with the Chinrest.
5. Ask the patient to:
 - Stand upright
 - Firmly grip the handle
 - Lean his/her chest lightly against the equipment
 - Position his/her feet slightly forward.
6. Adjust the chin support so that it rests snugly on the chin and secure the chin support firmly by turning the locking knob to the right.
7. Have the patient place his / her acanthion point on the chin support and tilt his / her head forward by 5 degree.

8. For **TMJ Open** imaging, have the patient:

- Open his/her mouth as wide as possible
- Keep his/her tongue positioned towards the palate
- Close his/her eyes.
- Be careful the chin not to touch any part of the equipment.

9. Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright. Ask the patient to remain still until scanning is complete.

CAUTION

To acquire the best image possible, ask the patient not to:




- Breathe or swallow saliva during image acquisition
- Move during image acquisition

Aligning the Laser Beams


Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.

- 1. Mid-sagittal plane laser beam:** Position the Mid-sagittal plane laser beam at the center of the patient's face to prevent magnifications of the left or right side in the final image.
- 2. Frankfurt plane laser beam:** Position the patient's head so that the Frankfurt plane is aligned with the Frankfurt plane laser beam.

- 3. Canine laser beam:** Have the patient smile to properly position the canine laser beam on the center of the patient's canine tooth.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head by clicking the temple support adjustment button.

2. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

3. Proceed to section **3.3 Initiating X-Ray Exposure**.

3.2.3 TMJ Close


Scan a **TMJ Open** image first and then scan a **TMJ Close** image. Guide the patient to close his/her lips for **TMJ Close** imaging while maintaining the same position used for **TMJ Open** imaging.

1. Shortly after **TMJ Open** imaging has been completed, you will receive the following message: "Do you want to take the **TMJ Close** image continuously?" Press the **OK** button to capture a **TMJ Close** image.

2. Adjust the chin support so that it rests snugly on the chin and secure the chin support firmly by turning the locking knob to the right.

3. Have the patient place his / her acanthion point on the chin support and tilt his / her head forward by 5 degree.

4. For **TMJ Close** imaging, Have the patient:


- Close his/her mouth
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes

5. The laser beams should be aligned in the same way as for **TMJ Open** imaging.

6. Proceed to section **3.3 Initiating X-Ray Exposure**.

3.2.4 Sinus


Positioning the Patient

1. Insert the Chin support into the Chinrest receptacle on the equipment.

Disinfect the Chinrest using an alcohol-based cleaning solution and wipe away all residues with a dry cloth before proceeding any further.

2. Open the **temple supports** by clicking the **temple support adjustment button**.

3. Guide the patient —facing the Chinrest — to the equipment. Have the patient stand upright in the center of the equipment.

4. Adjust the height of the equipment using the column up/down button or switch (optional) until the patient's chin rests comfortably on the Chinrest.


5. Ask the patient to:

- Stand upright
- Firmly grip the handle
- Lean his/her chest lightly against the equipment
- Position his/her feet slightly forward.

6. Adjust the chin support so that it rests snugly on the chin and secure the chin support firmly by turning the locking knob to the right.

7. Align the patient so that the area right under the lips touches the upper part of the chin support.

8. Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright.

9. Have the patient:

- Tilt his/her head back by 10-15°
- Close his/her mouth
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes

Ask the patient to remain still until scanning is completed.


To acquire the best image possible, ask the patient not to:

- **Breathe or swallow saliva during image acquisition**
- **Move during image acquisition**

Aligning the Laser Beams


Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.


1. Mid-sagittal plane laser beam: Position the Mid-sagittal plane laser beam at the center of the patient's face to prevent magnifications of the left or right side in the final image.

2. Frankfurt plane laser beam: The Frankfurt plane laser beam should go from the top of the ear to the tip of the nose when the patient's head is tilted back 10 – 15°.



Adjusting lever for the canine

3. Canine laser beam: Have the patient smile to properly position the canine laser beam on the center of the patient's canine tooth.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head by clicking the temple support adjustment button.

2. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

3. Proceed to section **3.3 Initiating X-Ray Exposure**.

3.3 Initiating X-Ray Exposure

The X-Ray exposure method and its sequences are the same for PANO Standard and PANO Special modes. The example used in this manual is the PANO Standard mode.

To emit the X-Ray, perform the following procedure.


- **Do not operate the PC during exposure. Failure to comply with this instruction may cause the system to malfunction.**
- **The operator shall observe the X-Ray safety regulations applicable to their area at all times during the operation of this equipment.**

If an emergency occurs during image acquisition, release the exposure switch to cease X-Ray emission.

1. Leave the X-Ray room and close the door. The operator must visually monitor the patient at all times during image acquisition.
2. Press and hold down the exposure switch until image acquisition is complete.

Exposure indicator
light turns orange
Orange: Exposure


During X-Ray irradiation, ensure that

- The LED lamp on top of the equipment turns orange to indicate X-Ray emission.
- The warning lamp outside the X-Ray room turns on.
- The sound (beep or music: optional) goes off.
- The radiation symbol on the upper left corner of GUI turns yellow to indicate X-Rays are being emitted.

3. The image appears in real time on the imaging GUI.

Ex) Magic PAN

Two images are obtained, with each in the **Magic PAN** and **Normal** respectively. You would click either button to compare them with each other.

The buttons:

Normal

Magic PAN

4. Click the **Save** button to save the image.

After Image Acquisition

After acquiring the image, complete the following steps:

- Loosen the temple supports to release the patient.
- Remove the hygiene cover from the bite block (for **standard PANO** mode only).
- Press **Return** to bring the Rotating Unit back to its initial position.

4 Acquiring CEPH Images

4.1 Setting the Exposure Parameters

NOTE

You can set the imaging parameters on either the touch screen or the imaging program running on the PC. They are synched in real time and display the same environmental settings.

Perform the following procedures to set the exposure parameters for the specific patient and capture mode.

Depending on the sensor type employed, one of three kinds of imaging S/W comes with the equipment for the CEPH mode examination.

Imaging Program

OP(1210 SGA sensor:
Oneshot Type)

Touch Screen

1. Click the **CEPH** button.
2. Select the scanning mode under **CEPH Examination**.

3. Select the FOV size for the specific imaging mode.

Model: OP / OS (One shot type sensor)

FOV (Inches)	Details	Model
12 x 10 (inches) 30.48x25.40 (cm)	Full size	OP
9 x 10 (inches) 22.86x25.40 (cm)	Uninterested region of the rear head is eliminated to minimize the X-Ray exposure area.	OP / OS
8 x 8 (inches) 20.32x20.32 (cm)	Child	OP / OS


4. EXPOSURE TIME

NOTE

Exposure time can be adjusted by resolution of 0.1 s in the range of 0.5 s to 1.2 s.

The Multi FOV is supported for the equipment with the one shot type sensor.


SC (Xmaru2301CF Sensor / Xmaru3001CF Sensor: Scan Type)

① Image Quality

Fast: scans faster with lower X-Ray dose than the **Normal** mode.

② CEPH Examination

Full Lateral (Optional): Scan a lateral view image fully.

5. The gender and age of the patient is automatically selected according to EasyDent / EzDent-i's patient information. Nevertheless, ensure that this information is correct.

Age Group		VATECH's Standard
Adult	Child	≤ 12
	Man	≥ 13
	Woman	

NOTE

A Child is defined as a person who is younger than 12 years old. If Child is selected, the image size and exposure dose are automatically reduced.

6. Select X-ray intensity.


NOTE

X-ray intensity (Hard, Normal, Soft) is according to the operator's decision.

Soft ≤ Normal ≤ Hard

7. A default value for the tube voltage (kVp) and current (mA) will be displayed based on the gender and X-ray intensity. If necessary, you may manually adjust these values using the arrows to the right of each value.

8. Click **CONFIRM** to accept the parameters.

When you click the **CONFIRM** button:

- The **READY** button will begin blinking to show that it is activate.
- The scan time and DAP (Dose Area Product) value are shown on the main display for the intended exposure.

Scan Time : 0.0
DAP : 0.00 uGy x m^2


9. Guide the patient to the equipment. Refer to section **4.2 Positioning the Patient**.

4.2 Positioning the Patient

Follow the procedure below to prepare and position the patient for CEPH imaging.

Before Positioning the Patient

- Ask the patient to remove all jewelry and metallic objects, such as earrings, hair pins, spectacles, dentures, and orthodontic appliances. These items may cause shadow images which could obscure diagnosis.
- It is strongly recommended that the patient wear a lead apron for protection against any possible scatter radiation.
- Adjust the column height to the patient's height using the up/down button or switch (optional).

NOTE

Correct positioning is an important factor in capturing the best possible image.

CAUTION

Ensure that the nasal positioner left unfolded, before adjusting the ear rods in the proper direction.

WARNING

Ensure that the patient is clear of all moving parts while adjusting the height of the equipment.


4.2.1 Lateral / Full Lateral

CEPH Lateral_ENG

*To link to the video-guided manual: scan the QR code with Smartphone or smart pad

Follow the procedure below to position the patient.


1. Widen the distance between the two ear rods.

NOTE

Porion position reference indicator enables the operator to easily identify the position of the porion on the image.

2. Guide the patient to the CEPH unit.


3. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).

4. Ask the patient to stand upright.

5. Make sure that the patient's shoulders are level and that their neck is relaxed.

6. The ear rods should fit properly along the patient's ear canals. The patient's Frankfurt Plane should be parallel to the floor.


7. Position the nasal positioner on the patient's Nasion point. The height of nasal positioner can be adjusted.

8. Ask the patient to swallow any saliva in his/her mouth and remain still until X-Ray exposure is complete.

After adjusting the height of the column to suit the patient, fit the ear rods along the patient's ear canals and adjust the nasal positioner.


9. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

10. Proceed to section 4.3 Initiating X-Ray Exposure.

4.2.2 Frontal (PA)

Follow the procedure below to correctly position the patient.

1. Turn the ear rods 90° clockwise from their initial position.


2. Widen the distance between the two ear rods.

3. The nasal positioner should be flipped to the side and upward to prevent it from obstructing the image acquisition.

Porion position reference indicator enables the operator to easily identify the position of the porion on the image.

4. Guide the patient to the CEPH unit.

5. Ask the patient to stand upright facing the sensor.


Make sure that the patient's shoulders are level and that his/her neck is relaxed.

6. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).


CAUTION

After adjusting the column to the patient's height, fit the ear rods along the patient's ear canals.

7. The ear rods should fit properly along the patient's ear canals. The patient's Frankfurt Plane should be parallel to the floor.

READY


8. Ask the patient to swallow any saliva in his/her mouth and remain still with his/her mouth closed until X-Ray exposure is complete.

9. Click the **READY** button after the patient has been positioned. No X-Ray will be emitted at this point.

10. Proceed to section **4.3 Initiating X-Ray Exposure**.


4.2.3 SMV

Follow the procedure below to prepare and position the patient.

1. Widen the distance between the two ear rods.

2. The nasal positioner should be flipped to the side and up to prevent it from obstructing the image acquisition.

3. Guide the patient to the CEPH unit.


4. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).

5. Ask the patient to stand upright facing the X-Ray tube.

Make sure that the patient's shoulders are level and that his/her neck is relaxed.

After adjusting the column to suit the height of the patient, fit the ear rods along the patient's ear canals.

6. Place the ear rods along the patient's ear canals. Make sure that the ear rods are comfortably, yet firmly in place.

7. Gently tilt the patient's head back until his/her Frankfurt plane is perpendicular to the floor, as shown below.


READY

8. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

9. Proceed to section **4.3 Initiating X-Ray Exposure**.

4.2.4 Waters View

Follow the procedure below to prepare and position the patient for Waters View mode image acquisition.

1. Widen the distance between the two ear rods.

2. Flip the nasal positioner to the side and up to prevent it from obstructing the image acquisition.

3. Guide the patient to the CEPH unit.

4. Ask the patient to stand upright facing the sensor. Make sure that the patient's shoulders are level and that his/her neck is relaxed.


5. Adjust the height of the unit to suit the patient by pressing the column up/down button or switch (optional).


CAUTION

After adjusting the column to suit the height of the patient, fit the ear rods along the patient's ear canals.

6. Place the ear rods along the patient's ear canals. Make sure that the ear rods are comfortably, yet firmly, in place.

7. Ask the patient to swallow any saliva in his/her mouth and tilt his/her neck back 30° - 40°, with the mouth closed, until X-Ray exposure is complete.


8. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.


9. Proceed to section **4.3 Initiating X-Ray Exposure**.

4.2.5 CARPUS

Attaching the CARPUS plate

1. Match the notch on the CARPUS plate to the end of nasal positioner.


2. Slide the CARPUS plate onto the nasal positioner.

3. Slide the locking bolt on the CARPUS plate towards the nasal positioner and turn it to secure the CARPUS plate in place.

4. Ensure that the CARPUS plate is firmly in place.

Positioning the Patient

1. Ask the patient to put his/her right hand flat on the CARPUS plate. It is important to ensure that the patient does not bend his/her fingers.

2. Have the patient close their eyes and remain still until scanning is complete.

Ask the patient not to block the ruler on the nasal positioner by his / her fingers. It could deteriorate image quality.

3. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

4. Proceed to section **4.3 Initiating X-Ray exposure**.

4.3 Initiating X-Ray Exposure

The X-Ray exposure method and process is the same for all CEPH modes. The example and images below are taken from an X-Ray performed in **CEPH Lateral** mode.

- Do not operate the PC during exposure. Failure to comply with this instruction may cause the system to malfunction.
- The operator shall observe the X-Ray safety regulations applicable to their area at all times during the operation of this equipment.


If an emergency occurs during image acquisition, release the exposure switch to cease X-Ray emission.

1. Leave the X-Ray room and close the door. The operator must visually monitor the patient at all times during image acquisition.

2. Press and hold down the exposure switch until image acquisition is complete.

Exposure indicator
light turns orange
Orange: Exposure

During X-Ray irradiation, ensure that


- The LED lamp on top of the equipment turns orange to indicate X-Ray emission.
- The warning lamp outside the X-Ray room turns on.
- The sound (beep or music: optional) goes off.
- The radiation symbol on the upper left corner of GUI turns yellow to indicate X-Rays are being emitted.

The image appears in real time on the imaging GUI.

3. Click the **Save** button to save the acquired image. If automatic save has been selected as a default setting, the image will be saved automatically.

Sample images: with One shot type sensor

Lateral (Unit: inches)

**Full Lateral : FOV 12 x 10 (inches)
30.48x25.40 (cm)**

**Lateral : FOV 9 x 10 (inches)
22.86x25.40 (cm)**

**FOV 8 x 8 (inches)
20.32x20.32 (cm)**

Carpus

After Image Acquisition

After the image has been acquired, perform the following tasks:

- Fold away the nasal positioner.
- Loosen the ear rod supports and remove them from the patient's ears.
- Release the patient.

5 Acquiring CT Images


5.1 Setting the Exposure Parameters

NOTE



You can set the imaging parameters on either the touch screen or the imaging program running on the PC. They are synched in real time and display the same environmental settings.

Perform the following procedures to set capture parameters for the specific patient and capture mode.

Ex: with Xmaru1524CF
Master Plus sensor

Touch Screen

NOTE


The FOV size is determined by the CT sensor attached on the equipment.

1. Select the FOV size.

2. Select the ROI to be captured under VERTICAL POSITION and HORIZONTAL POSITION. Those information selected will be displayed on the guidance image.

For the FOV 50 X 50, you may select a specific tooth position to be captured on the full arch guidance image.

Guidance Image	FOV (mm)	ROI
	160 X 100	Mn - Center
	150 X 150	Mn - Center
	120 X 90	Mn - Center
	100 X 80	Mn - Center
	80 X 80	Mn - Center
	80 X 50	Mn - Center
	50 X 50	Mn / Left Central Incisor

3. Select the image quality mode.

4. Select the voxel resolution.

5. If **Metal Artifact Reduction** is necessary, click **Apply** under **Metal Artifact Reduction**. Metal Artifact Reduction reduces the appearance of metal in imaging. This function is most effective when there are fewer than 3 metal artifacts.

NOTE

The MAR (Metal Artifact Reduction) function doubles image reconstruction time.

6. Gender selection will be made automatically according to the patient information registered in EasyDent / EzDent-i. Nevertheless, ensure that this information is correct.

Age Group		VATECH's Standard
Child		≤ 12
Adult	Man	≥ 13
	Woman	

7. Select X-ray intensity.


NOTE

X-ray intensity (Hard, Normal, Soft) is according to the operator's decision.

Soft \leq Normal \leq Hard

8. A default value for the tube voltage (kVp) and current (mA) will be displayed based on the gender and X-ray intensity. If necessary, you can perform manual adjustments by clicking the arrows to the right of each number.

9. Click **CONFIRM** to accept these parameters.

NOTE

Please wait for a moment, while the rotating unit moves to its initial scanning position.

When you click **CONFIRM** button:

- The **CONFIRM** button turns to the **SCOUT** and then the **READY** button is enabled and flashes. This means that the equipment is ready for X-Ray exposure.

NOTE

SCOUT capture is featured in FOV size 50 X 50 and 80 X 50. In case of abnormal patient's arch, the chinrest height adjustment through the SCOUT feature enables the user to acquire better image.

When the SCOUT image needs to be acquired, position the patient first, followed by clicking the SCOUT button.

For details about imaging in the SCOUT mode, refer to the section "5.3 Acquiring a Scout image".

- The Rotating unit will move to its initial scanning position.
- The Chinrest will adjust itself according to the mode selected.
- Three laser beams (**CT horizontal, Mid-sagittal and CT vertical**) will activate to facilitate patient alignment. These beams will disappear after 20 minutes or if the **READY** button is clicked.

NOTE

Laser beam On / Off buttons

: Imaging Program

: Control Panel

Scan Time : 0.0
DAP : 0.00 uGy x m²

- The scan time and DAP (Dose Area Product) value are shown on the main display for the intended exposure.


10. Guide the patient to the equipment. Position the patient. Refer to section **5.2 Positioning the Patient**.

5.2 Positioning the Patient

Follow the procedure outlined below to properly position the patient for CT scanning.

Before Positioning the Patient

- Ask the patient to remove all jewelry and metallic objects, such as earrings, hair pins, spectacles, dentures, and orthodontic appliances. These items can cause shadow images which could obstruct diagnosis.
- It is strongly recommended that the patient wear a lead apron to protect against scatter radiation.
- Adjust the column height to the patient's height using the up/down button or switch (optional).


Correct positioning is an important factor in capturing the best possible image. Proper positioning will reduce the appearance of the cervical spine in the image.

Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.

Positioning the Patient

1. Insert the bite block into the Chinrest receptacle on the equipment.


2. Place a new hygienic cover over the bite block.

The hygiene cover for the bite block is intended for single use only. Always replace the hygiene cover for each new patient.

Disinfect the Chinrest and bite block using an alcohol-based cleaning solution and wipe away all residue with dry cloth before each new patient.

3. Open the temple supports by clicking the Temple Support Wheel Adjustment button.

4. Guide the patient —facing the Chinrest — to the equipment.


5. Adjust the height of the column using the column up/down button or switch (optional) until the patient's chin rests comfortably on the Chinrest.

6. Ask the patient to:

- Stand upright
- Firmly grip the handles on both sides of the unit
- Lean his/her chest lightly against the equipment
- Position his/her feet slightly forward and open stance

Ensure that the patient's shoulders remain level and that his/her neck is relaxed. The cervical spine should be straight and upright.

7. Have the patient:

- Bite the bite block along its grooves using his/her upper and lower incisors
- Close his/her lips
- Keep his/her tongue pressed against his/her palate
- Close his/her eyes

The patient's occlusion line should be parallel to the floor.

Ask the patient to remain still until scanning is complete.


To acquire the best image possible, ask the patient not to:

- Breathe or swallow saliva during image acquisition
- Move during image acquisition

Aligning the Laser Beams


Be careful not to project the laser beams directly into the patient's eyes as this could severely damage the patient's vision.

The 3D CT image centers the FOV on the intersection of the three laser beams (CT horizontal, Mid-sagittal and CT vertical).


1. CT horizontal laser beam: Position the CT horizontal laser beam at the center of the FOV area.

2. Mid-sagittal plane laser beam: Position the Mid-sagittal plane laser beam at the center of the FOV area.

3. CT vertical laser beam: Position the CT vertical laser beam at the center of FOV from the patient's lateral.

Completing Patient Positioning

1. After confirming patient positioning and beam alignment, adjust the temple supports to fit snugly on both sides of the patient's head using the temple support adjustment button.

2. Click the **READY** button after the patient has been properly positioned. No X-Ray will be emitted at this point.

If, for the FOV size of 50 x 50 or 80 X 50, the **SCOUT** image is to be acquired, click the **SCOUT** button.

NOTE

SCOUT capture is featured in FOV size 50 X 50 and 80 X 50.

In case of abnormal patient's arch, the chinrest height adjustment through the SCOUT feature enables the user to acquire better image.

For details about imaging in the SCOUT mode, refer to the section "5.3 Acquiring a Scout Image".


3. Proceed to the section **5.4 Initiating X-Ray Exposure**. But if the **SCOUT** feature is selected, continue to the next section **5.3 Acquiring a SCOUT Image**, to the section **5.4 Initiating X-Ray Exposure**.

5.3 Acquiring a SCOUT Image

In general, the imaging trajectory is predefined along the normal arch shape. In case of abnormal patient's arch, the chinrest height adjustment through the SCOUT feature enables the user to acquire better image.


5.3.1 Starting imaging with the SCOUT feature

SCOUT Imaging Flow Chart

Imaging and Coordinate compensation

The parameter setting for exposure and patient positioning are the same as those described in the section 5.1 and 5.2. For the SCOUT feature-specific imaging, do take the steps in the following manner.

1. When the parameter setting and patient position are done, click the **SCOUT** button.
2. Press and hold down the exposure switch to acquire the **SCOUT** image.

3. Guide the patient out of the equipment.

Always perform the chinrest coordinate adjustment after the patient is cleared from the equipment.

Failure to do so could result in the physical injury, due to the vertical movement of chinrest to new position reflected after the 'CONFIRM' button on the SCOUT Viewer is clicked.

4. Once the imaging has been completed, the resulting image with the SCOUT feature appears.

The green-colored horizontal line appears on the **SCOUT** image, when the left button of the mouse is clicked at the point which you want to move to the center of the image.

To adjust the image contrast for making the compensating point distinct, you may adjust the black and white level with the Adjust control on the right of image.

For details about the SCOUT Viewer, refer to the section "5.3.2 SCOUT Viewer".

5. Upon clicking the **CONFIRM** button, the **SCOUT** view screen closed, and the chinrest moves to newly compensated position vertically.

6. Guide the patient into the equipment again and position the patient.


7. Click the **READY** button from the GUI.

8. Press and hold down the exposure switch to acquire the CT image: refer to the section **5.4 Initiating X-Ray Exposure**.

If the acquired image is unsatisfactory, you may repeat the steps above with a different coordinate.

5.3.2 SCOUT Viewer

1. Raw View Window

- Displays the acquired image and the guideline.
- The green-colored horizontal line appears when the mouse clicked at the point to be moved to on the image.

2. Adjust

- Adjust the current contrast level

White: make the current white level whiter

Black: make the current dark level darker

3. Initialize

- When clicked, returns to the initial state of the completion of the **SCOUT** imaging.

4. CONFIRM

- When clicked, the **SCOUT** view screen closed and the chinrest moves to newly compensated position vertically.


5.4 Initiating X-Ray Exposure

- **Do not operate the PC during exposure. Failure to comply with this instruction may cause the system to malfunction.**
- **The operator shall observe the X-Ray safety regulations applicable to their area at all times during the operation of this equipment.**

1. Leave the X-Ray room and close the door. The operator must visually monitor the patient at all times during image acquisition.
2. Press and hold down the exposure switch until image acquisition is complete.

Exposure indicator
light turns orange
Orange: Exposure

During X-Ray irradiation, ensure that

- The LED lamp on top of the equipment turns orange to indicate X-Ray emission.
- The warning lamp outside the X-Ray room turns on.
- The sound (beep or music: optional) goes off.
- The radiation symbol on the upper left corner of GUI turns yellow to indicate X-Rays are being emitted.

The image appears in real time on the imaging GUI.

If an emergency occurs during image acquisition, release the exposure switch to cease X-Ray emission.

3. The captured image is automatically reconstructed and converted into a DICOM file. Click the Save button to store the file in the patient database (If automatic save is selected in the program settings, the image will be saved automatically).

Once acquisition is complete, the image is automatically transferred to **EasyDent / EzDent-i**.

For more details about 2D or 3D viewer, refer to EasyDent / EzDent-i, and Ez3D plus / Ez3D-i user manual.

After Image Acquisition

Perform the following tasks after acquiring a CT image:

- Loosen the temple supports and release the patient.
- Remove the hygiene barrier from the bite block.
- Press Return to bring the Rotating Unit back to its initial position.

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Copyright © 2014 par VATECH
Copyright © 2014 por VATECH